

ResMed

Actual Size:
148 mm (Width)
x 110 mm (Height)
or
5^{3/4}" (Width)
x 4^{1/4}" (Height)




Ruler is for scaling purposes.


How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size

AirTouch™ N30i

Nasal mask






ResMed



AirTouch™ N30i
Nasal mask

AirFit™ N30i
Nasal mask

AirFit™ P30i
Nasal pillows mask

Frame Fitting Instructions



 Small Frame	 Standard Frame
0-150mm	>150mm
0-6"	>6"